


# HIOKI Powder Impedance Measurement System

## All-in-One Solution for Solid-State Batteries and Dry Processes




### Simultaneously measure Impedance, Thickness, and Pressure


- Accurate Analysis**  
 Calculate **bulk density** and **ionic conductivity** with precision
- All-In-One Glove Box Operation**  
 From loading to pressing and measurement, all tasks completed safely inside.
- Enhanced Safety**  
 Prevents hydrogen sulfide gas leakage and preserves material integrity.
- Time-Saving Efficiency**  
 No sample removal required: Streamline the entire process
- Optimized Testing**  
 Measure multiple conditions on a single sample continuously
- Wide frequency range**  
 4 Hz to 5 MHz for better analysis
- Automatic Calculations**  
 Easily calculate total resistance (R) with Nyquist plot fitting.
- Outputs Electrical & Physical Parameters**  
 Including conductivity, resistivity, bulk density, packing density, and porosity.
- Graphing Capabilities**  
 Visualize results with built-in graphing functions for in-depth analysis.



### System Configuration



### Real data: Measuring Differences Between Manufacturers and Grain Sizes For Sulfide-Based Solid Electrolytes



- Measurement samples**
  - Li<sub>6</sub>PS<sub>4</sub>Cl argyrodite
  - Manufacturer (particle size): Company A (1 μm), Company B (8 μm)
  - Mass: 200 mg
- Measurement conditions**
  - Measurement frequency: 4 Hz to 5 MHz
  - Signal level: 0.1 V (constant voltage)
  - Electrode diameter: φ10 mm

- When the press pressure increases, the overall resistance R decreases.
- Differences between sample manufacturers and grain sizes can be identified.

