Evaluating the Electron Conductivity of an LIB Electrode Slurry

Introduction of the Slurry Analytical System to estimate the dispersion condition by evaluating the electronic conductivity of a slurry.

Applications

Searching for materials, Optimizing manufacturing conditions, Screening slurries based on performance

Background

Battery electrode materials are used in the form of slurries that are diffused in a dispersion medium. Since the state of dispersion in a slurry is unstable, it's important to evaluate dispersion conditions as part of the quality control process. In addition, since the slurry manufacturing process occurs at the beginning of battery manufacturing, screening needs to be carried out at an earlier stage.

This Application Note introduces an application in which dispersion conditions are estimated by evaluating electron conductivity in an LIB electrode slurry.

For more information about the principles underlying the Slurry Analytical System, please see this FAQ. (https://www.hioki.com/global/support/faq/detail/id_n1265055)

1. Differences in mixing time

Objective

To optimize mixing time

Slurry production conditions

Solid content ratio	Active material: 95 % Conductive auxiliary agent: 2.5 % Binder: 2.5 %	
Mixing process	Mix for 0, 1, 3, or 5 min. with a solid content concentration of 70%	
Viscosity adjustment process	Dilute with NMP until the solid content concentration reaches 50%	

Analysis results



Uniformity
0.82
0.8
0.78
0.76
0.76





© 2023 HIOKI E.E. CORPORATION application_slurry_case_E2_230203

All information correct as of February 3, 2023. Contents are subject to change without notice. Note: Company names and product names appearing in this brochure are trademarks or registered trademarks of various companies.

Observations

With increasing mixing time, the Rratio decreased, and the DCR increased. These results suggest that shorter mixing time is associated with conductive material network development, which causes the slurry's overall resistance to decrease. By contrast, uniformity increased with mixing time. This result suggests that longer mixing time is favorable for slurry dispersibility. When these results are combined, we can surmise that mixing time induces the following states in the conductive materials in the slurry, and that the mixing time of 3 min or longer yields a comparatively favorable slurry.



Mixing time of 0 to 1 min Conductive materials are in a state of agglomeration, resulting in short, and thick conductive paths.



Mixing time of 3 to 5 min Conductive materials are in a state of dispersion, resulting in long, and thin conductive paths.

2. Differences in conductive auxiliary agent quantity

Objective

To verify the effect of adding conductive auxiliary agent as a way to reduce electrode resistance To verify the relationship between slurry impedance and electrode resistanc

Slurry production conditions

Solid content ratio	Active material: 18 g, Conductive auxiliary agent: 0.25 g, 0.5 g, 0.75 g, Binder: 0.5 g
Mixing process	Mix for 9 min with a solid content concentration of 70%
Viscosity adjustment process	Dilute with NMP until the solid content concentration reaches 50%

Analysis results



ΗΙΟΚΙ

© 2023 HIOKI E.E. CORPORATION application_slurry_case_E2_230203

All information correct as of February 3, 2023. Contents are subject to change without notice.

Note: Company names and product names appearing in this brochure are trademarks or registered trademarks of various companies.

Observations

The Rratio rose, and the DCR fell, when the quantity of conductive auxiliary agent was increased. These changes, which suggest that adding conductive auxiliary agent causes conductive material networks to develop in the slurry, lowering its overall resistance to decrease, are in line with the theory underlying slurry analysis. By contrast, uniformity with the 0.75 g of auxiliary agent was lower than with the 0.5 g. We can surmise that conductive material networks formed more readily with the 0.75 g of auxiliary agent, but that dispersion was insufficient.

We formed electrode sheets with these slurries and analyzed them with the Electrode Resistance Measurement System RM2610. The results indicated that both the composite layer resistance and the interface resistance fell as the quantity of conductive auxiliary agent increased. These results suggest that the slurry's electron conductivity has been inherited by the electrode sheets' resistance characteristics, indicating a good application process.

3. Relationship between electron conductivity and dynamic viscoelasticity resulting from differences in dispersant quantity

provided by Shin-Etsu Chemical Co., Ltd.

Objective

To verify the effect of dispersant quantity on a slurry using both in terms of impedance and rheology

Slurry production conditions

Solid content ratio	Active material: 96 %, Conductive auxiliary agent: 2 %, Binder: 2 %, Dispersant: 0.08% to 0.16% (4% to 8% of the conductive auxiliary agent quantity)
Conductive auxiliary agent slurry preparation	Preparation of a conductive auxiliary agent slurry containing only electrode auxiliary agent, dispersant, and NMP
Electrode slurry preparation	Preparation of an electrode slurry by adding active material and binder

Analysis results



Observations

Adding dispersant from 4% to 5% of the conductive auxiliary agent quantity

Rratio increased abruptly, and the difference between G' and G' decreased.

These results suggest that adding dispersant eliminates agglomeration of the conductive auxiliary agent and that the conductive auxiliary agent networks formed are functioning as conductive paths.

Adding dispersant from 5% to 8% of the conductive auxiliary agent quantity

We can surmise that adding even more dispersant will increase fluidity and improve the coating properties due to the slurry's favorable G" value, while the reduction in Rratio means that the conductive auxiliary agent networks will become even thinner, leading to breakages.



© 2023 HIOKI E.E. CORPORATION application_slurry_case_E2_230203

Services and equipment used

SLURRY ANALYTICAL SYSTEM		HIOKI Service (license purchase required for use)
LCR METER	IM3536	HIOKI Product
ANALYTICAL SOFTWARE	SA2632	HIOKI Free Software
MEASUREMENT SOFTWARE	SA2633	HIOKI Free Software
ELECTRODE CELL	SA9001	HIOKI Product
TEST FIXTURE	SA9002	HIOKI Product
ELECTRODE RESISTANCE MEASUREMENT SYSTEM	RM2610	HIOKI Product

Sein,



© 2023 HIOKI E.E. CORPORATION application_slurry_case_E2_230203

All information correct as of February 3, 2023. Contents are subject to change without notice. Note: Company names and product names appearing in this brochure are trademarks or registered trademarks of various companies.